Vangerowstraße 14
69115 Heidelberg
+49 (0)6221 - 90 50 50

www.finescience.de
europa@finescience.de

DEUTSCH	(ab S. 2)
ENGLISH	(from p. 4)
FRANÇAIS	(à partir de la p. 6)
ITALIANO	(da p. 8)
ESPAÑOL	(desde p. 10)

FST # 24201-00

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
Laboratory Animal Microtattoo System

The Laboratory Animal Microtattoo System is a quick, efficient and affordable method of permanently identifying laboratory animals, particularly mice and rats as early as one day old.

Designed and tested by laboratory animal technicians, the system utilizes a hypodermic needle at one end of a stainless steel forceps and a small metal cup containing green tattoo paste on the opposite forceps. By first squeezing the forceps together to pick up tattoo paste on the needle, and then squeezing the forceps so that the needle penetrates the toe pad, foot pad, ear or tail of the animal, a small identifying mark is created.

A convenient magnifier is attached to the instrument to assist with injecting the tattoo paste clearly and accurately.

Items required:

- Microtattoo instrument ⁺
 - Green tattoo paste ⁺
 - Sterile hypodermic needle Recommended size:
 - Adult mice: 27 ga x 1/2" ⁺
 - Young mice: 30 ga x 1/2"
 - Rats and other small animals: 25 ga x 5/8"
 - Aluminum ink planchette ⁺
 - Alcohol (used for sterilization)
- ⁺ *Included with FST system.*

Ink may be sterilized by autoclave.

Before starting, determine the grid location pattern to be used, and then follow the numbered directions. A sample grid pattern is provided opposite.

1. Place a small amount of tattoo paste on the aluminum planchette.
2. Sterilize the needle with alcohol.
3. Dip the reservoir end of the instrument in the tattoo paste and fill the reservoir. *
4. Squeeze the instrument ends together so that the needle goes into the reservoir and acquires some paste.
5. Immobilize the animal according to your normal method.
6. Squeeze the instrument ends together so that the needle completely penetrates the toe pad, foot pad, ear or tail of the animal. Tattoo paste will be deposited in the area penetrated by the needle.
7. After use, the instrument should be thoroughly cleaned with warm water and soap.


* Alternate method to fill reservoir:

Instead of directly filling the reservoir, you can use a modified pipette tip.

Cut the smaller end of a standard pipette tip to snugly fit into the reservoir. Some experimentation will be required to achieve the optimum fit and to make sure the needle does not catch once inserted.

Fill the modified pipette tip with tattoo paste directly from the tube or with a syringe and screw it into the reservoir opening.

You can fill more than one pipette tip in advance, and replace them as needed, avoiding the need to clean the tool constantly when working with large numbers of animals.

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